

Magnet Sales Information


300 MHz / 154 mm

UltraShield™ Plus – Super Wide bore - Long hold time



Magnet System 300/154 US PLUS™ LH
Part Number Z106379

Geometrical Dimensions



Z1044563

Geometrical Dimensions

| | | Description |
|--------|---------|--|
| C = | 3410 mm | Operational ceiling height |
| D-RT = | 154 mm | Diameter room temperature bore tube |
| D1 = | 800 mm | Diameter cryostat upper part |
| D2 = | 850 mm | Diameter cryostat bottom plate |
| H1 = | 1140 mm | Height of cryostat from bottom flange – upper flange |
| H2 = | 1526 mm | Height of cryostat from bottom flange to Helium tower Minimum height for transportation |
| H3 = | 1699 mm | Height of cryostat from bottom flange to Helium manifold |
| S = | 1050 mm | Height between floor and magnet bottom flange |


System Data

| | |
|---|-----------------------|
| Minimum operational ceiling height (Helium transfer line 29085) | 3130 mm |
| Minimum ceiling height with standard Helium transfer line 53962 | 3410 mm |
| Required space (footprint, width x depth) | ~ 1.49 m ² |
| System weight (empty, without magnet stand) | 534 kg |
| Magnet stand | 115 kg |
| System weight (filled completely, with magnet stand ADI/EMI) | 780/765 kg |

NMR Magnet Specifications

| | |
|--|--|
| Type | BZH 300'154 US PLUS™ |
| NMR-frequency (¹ H) | 300 MHz |
| Operating field | 7.05 Tesla |
| Field stability (guaranteed value in persistent mode) | < 20 ppb/hr (< 6 Hz/hr) |
| Axial range with homogeneity better than 10 ppm | ~ 60 mm |
| Radial stray field (horizontal distance of the 0.5 mT (5G line from the magnetic centre) | < 1.00 m |
| Axial stray field (vertical distance of the 0.5 mT (5G line from the magnetic centre) | < 2.00 m |
| Cryo shims | X, Y, Z, Z ² , Z ³ , XZ, YZ, X-Y, X ² -Y ² |

Fringe Field Plot



| Fringe field contour | Radial [R] | Axial [Z] |
|---------------------------------|-----------------|-----------------|
| 200 mT (Directive 2004/40/EC) | Inside cryostat | Inside cryostat |
| 5 mT | 0.65 m | 1.18 m |
| 3 mT | 0.72 m | 1.32 m |
| 1 mT | 0.88 m | 1.69 m |
| 0.5 mT (5 Gauss) | 1.00 m | 2.00 m |
| 0.2 mT | 1.29 m | 2.53 m |
| 0.1 mT | 1.63 m | 3.06 m |
| 0.05 mT (~Earth magnetic field) | 2.10 m | 3.74 m |

Cryostat Specifications

| | |
|---|-------------------------------|
| Type | D 345/154 US PLUS™ |
| Room temperature bore | 154 mm |
| Approx. Helium evaporation rate under stabilized conditions (T=20°C, p=1030 mbar) | < 28 ml liquid Helium/hour |
| Liquid Helium refill volume/total volume | ~ 60/101 litres |
| Helium hold time | > 90 days |
| Approx. Nitrogen evaporation rate under stabilized conditions (T=20°C, p=1030 mbar) | < 300 ml liquid Nitrogen/hour |
| Liquid Nitrogen refill volume/total volume | ~ 101/126 litres |
| Nitrogen hold time | > 14 days |


Accessories

| | | |
|---|----------|------------|
| Magnet stand F85-1050 EMI (height 1050 mm) E lastomer I solators with vertical damping Frequencies damped > 14 Hz / resonance frequency = 9.5 Hz | Standard | Z109694 |
| Magnet stand F85-1050 ADI (height 1050 mm) A ir Spring and D amped I solators with vertical damping Frequencies damped > 3.8 Hz / resonance frequency = 2.6 Hz | Optional | Z109426 |
| Magnet stand F85-1050 API (height 1050 mm) A ir P iston and Damped I solators with vertical and horizontal damping Frequencies damped > 3.8 Hz / resonance frequency = 2.6 Hz | Optional | On request |
| Nitrogen level sensor for BSMSII (SCB3) (not working with BSNL) | Optional | Z122400 |

Equipment for Cryogen Transfer

| | | |
|---|--------|-------|
| Helium transfer line* D3xx (1455/2060/655) | AH0070 | 53962 |
| Helium transfer line* with bendable extensions (1455/2060/380) for minimum operational ceiling height (3000 mm) | | 29085 |

Top View



Geometrical Dimensions

| | | |
|---------------------------------------|----|---------|
| Width of magnet stand | W1 | 1285 mm |
| Depth of magnet stand | W2 | 1157 mm |
| Diameter of magnet stand = 2 x radius | D3 | 1430 mm |

Transport

| | | |
|--|-----------|--------------------------------|
| Overall system dimensions for transportation | | |
| Magnet system box | L x D x H | 105x 126 x 180 cm ³ |
| Magnet stand box | L x D x H | 167 x 88 x 101 cm ³ |
| Minimum system dimensions of magnet, unpacked (without manifold) | | Ø 85 cm, H 153 cm |
| System weight for transportation | | ~ 710 kg |
| Weight magnet stand box (ADI/EMI) | | ~ 198/182 kg |

Installation

| | |
|--|-------------------------|
| Liquid Nitrogen needed for installation (cooling down) | 480 litres |
| Liquid Helium needed for installation (cooling down, energizing, shimming, filled up completely) | 200 litres |
| Liquid Helium needed after a quench | 110 litres |
| Nitrogen gas for flushing | 1 cylinder 50 l/200 bar |
| Helium gas for flushing | 1 cylinder 50 l/200 bar |

* A detailed description of the marked objects can be found in chapter „Accessories“.



Bruker BioSpin

- **Bruker BioSpin**
your solution partner

Bruker BioSpin provides a world class, market-leading range of analysis solutions for your life and materials science needs

- **Bruker BioSpin Group**

sales@bruker-biospin.com
www.bruker-biospin.com